

**PurePRP® SupraPhysiologic Concentrating System
GenesisCS Component Concentrating System**

**Date: May 2025
Instruction for use**

ATTENTION OPERATING SURGEON

NOTE: DEVICE IS FOR SINGLE USE ONLY. Discard the entire disposable system after one use, using an acceptable disposal method for products potentially contaminated with blood.

DESCRIPTION

1. The PurePRP® Supraphysiologic Concentrating System is manufactured by EmCyte Corporation. The kit prepares platelet rich plasma from a small sample of blood at the point of care. The system contains syringes, needles and the concentrating device accessories.

MATERIALS

2. The materials used are syringes, needles, tubing, connectors, and concentrating devices. The materials consist of medical grade polymers, elastomers and stainless steel that are suitable for use in medical devices. All components in this system are packaged, labeled and sterilized as indicated by the manufacturer's labeling. All components in this system are latex-free.

INDICATIONS FOR USE STATEMENTS

3. The PurePRP® Supraphysiologic Concentrating System is designed to be used for the safe and rapid preparation of autologous platelet rich plasma (PRP) from a small sample of blood at the patient's point of care. The PRP is mixed with autograft and allograft bone prior to application to an orthopedic site to improve bone graft handling characteristics.
4. The safety and effectiveness of this device for in vivo indications for use, such as bone healing and hemostasis, have not been established.
5. The PRP prepared by this device has not been evaluated for any clinical indications.
6. The PRP prepared by this device is NOT indicated for delivery to the patient's circulatory system.

USER POPULATION

7. The intended user population is medical professionals who are licensed or certified in clinical practice. The operational context of the device requires users to be trained on aseptic technique and understand blood components. The surgeon is to be thoroughly familiar with the equipment and the surgical procedure prior to using this device.

DEVICE USE ENVIRONMENT

8. The device is intended to be used in a health care setting such as a surgery room, clinic or outpatient care center.

WARNING AND PRECAUTIONS

9. Use proper safety precautions to guard against needle sticks.
10. Follow manufacturer instructions when using centrifuge. Use only EmCyte provided general purpose centrifuge. Outcomes using centrifuges from other manufacturers are unknown.
11. Do not use sterile components of this system if package is opened or damaged.
12. Single use device. Do not reuse. Do not attempt to clean or re-sterilize this product.
13. Do not use after expiration date.
14. Use prepared PRP within 4 hours after drawing blood according to current AABB guidelines.

POSSIBLE RISKS

15. The patient is to be made aware of the general risks associated with whole blood aspiration. These risks include, but are not limited to: hemorrhage, seroma formation, infection, and/or persistent pain at the site of aspiration.
16. Reuse may be a potential biohazard

POSSIBLE ADVERSE EFFECTS

17. Damage to blood vessels, hematoma, delayed wound healing and/or infection is associated with blood draw, and/or surgical procedure.
18. Temporary or permanent nerve damage that may result in pain or numbness is associated with blood draw, and/or surgical procedure.
19. Early or late postoperative infection is associated with surgical procedure.
20. Pain associated with site of whole blood harvest.

STERILITY

21. The PurePRP® SupraPhysiologic Concentrating System kits are sterilized by ETO exposure. Do not use any component from an opened or damaged package. Do not re-sterilize. Discard if kit packaging is damaged or open.

INSTRUCTIONS FOR USE FOR 60mL SYSTEM

PREPARATION PROTOCOL:

22. NOTE: Use standard sterile aseptic technique throughout the following procedure. Always swab needle-less ports with alcohol before and after accessing.
23. WHOLE BLOOD DRAW: Attach the sterile filter needle onto the sterile 60mL syringe. Draw 8mL of Citrate Anticoagulant into the 60mL syringe. Remove the filter needle from the syringe. Attach the butterfly needle onto 60mL syringe and prime the needle with the anticoagulant. Slowly draw 52mL of whole blood from the patient filling the syringe to 60mL. Gently, but thoroughly mix the blood and anticoagulant upon collection to prevent coagulation.

CONCENTRATING PROTOCOL:

24. LOAD: IMPORTANT: Attach sterile non-vented clear cap to the bottom port of the device. **STERILE CLEAR CAP MUST BE ATTACHED TO THE BOTTOM PORT BEFORE CENTRIFUGATION.** Slowly add the anticoagulated whole blood through the top port of the Concentrating Device.
25. BALANCE: Make sure the counterbalance device contains the same amount of volume as the Concentrating Device. Then place them directly opposite to each other in the centrifuge rotor buckets. Close the lid.
26. FIRST SPIN:
 - a. Sapphire Series Centrifuge: PUREPRP 60 SPIN 1.
 - b. Platinum Series Centrifuge: PUREPRP SP SPIN 1.
 - c. Executive Series Centrifuge: 2.0 minutes and 4400 RPM.
 - d. Press the start button. Once the centrifuge stops, remove the **Device**.
27. FIRST EXTRACTION & TRANSFER: Attach the sterile 60mL syringe to the top port.
 - a. **LP-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate the remaining plasma buffycoat, leaving the RBC layer behind.
 - b. **MR-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate collecting 0.5mL of RBC buffycoat.
 - c. **NR-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate collecting 1mL of RBC buffycoat.
28. SECOND SPIN: Place Concentrating Device back into the centrifuge rotor bucket directly opposite the counterbalance device. Close the lid.
 - a. Sapphire Series Centrifuge: PUREPRP 60 SPIN 2.
 - b. Platinum Series Centrifuge: PUREPRP SP SPIN 2.
 - c. Executive Series Centrifuge: 4.0 minutes and 4400 RPM.
 - d. Press the start button. Once the centrifuge stops, remove the **Device**
29. SECOND EXTRACTION: Remove the clear cap from the bottom port. Using the 60mL syringe, aspirate plasma from the bottom port leaving 7mL or the desired amount in the device.
30. RESUSPEND THE PRP: Gently swirl the Concentrating Device to re-suspend the platelet concentrate into the plasma.
31. EXTRACT PRP: Attach a sterile 12mL syringe to the bottom port and tilt to immerse the aspirating pipe, then aspirate the platelet rich plasma through the open port of the aspirating pipe. Remove sterile syringe and apply a sterile cap.

**LP – LEUKOCYTE POOR
MR – MONOCYTE RICH
NR – NEUTROPHIL RICH**

INSTRUCTIONS FOR USE FOR 30mL SYSTEM

PREPARATION PROTOCOL:

32. NOTE: Use standard sterile aseptic technique throughout the following procedure. Always swab needle-less ports with alcohol before and after accessing.
33. WHOLE BLOOD DRAW: Attach the sterile filter needle onto the sterile 30mL syringe. Draw 3mL of Citrate Anticoagulant into the 30mL syringe. Remove the filter needle from the syringe. Attach the butterfly needle onto 60mL syringe and prime the needle with the anticoagulant. Slowly draw 27mL of whole blood from the patient filling the syringe to 30mL. Gently, but thoroughly mix the blood and anticoagulant upon collection to prevent coagulation.

CONCENTRATING PROTOCOL:

34. LOAD: IMPORTANT: Attach sterile non-vented clear cap to the bottom port of the device. **STERILE CLEAR CAP MUST BE ATTACHED TO THE BOTTOM PORT BEFORE CENTRIFUGATION. WITH TOP VENT OPEN**, slowly add the anticoagulated whole blood through the top port of the Concentrating Device. **THEN CLOSE VENT.**
35. BALANCE: Make sure the counterbalance device contains the same amount of volume as the Concentrating Device. Then place them directly opposite to each other in the centrifuge rotor buckets. Close the lid.
36. FIRST SPIN:
 - a. Sapphire Series Centrifuge: PUREPRP 30 SPIN 1.
 - b. Platinum Series Centrifuge: PUREPRP SP SPIN 1.
 - c. Executive Series Centrifuge: 1.0 minutes and 4400 RPM.
 - d. Press the start button. Once the centrifuge stops, remove the **Device**.
37. FIRST EXTRACTION & TRANSFER: Attach the sterile 30mL syringe to the top port.
 - a. **LP-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 30mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate the remaining plasma buffycoat, leaving the RBC layer behind.
 - b. **MR-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 30mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate collecting 0.5mL of RBC buffycoat.
 - c. **NR-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate collecting 1mL of RBC buffycoat.
 - d. **Remove CLEAR CAP** from bottom port and inject the 3mL syringe buffycoat solution, followed by the platelet plasma suspension (PPS). **Then SNUGGLY PLACE STERILE CLEAR CAP ON THE BOTTOM PORT!!!**
38. SECOND SPIN: Place Concentrating Device back into the centrifuge rotor bucket directly opposite the counterbalance device. Close the lid.
 - a. Sapphire Series Centrifuge: PUREPRP 30 SPIN 2.
 - b. Platinum Series Centrifuge: PUREPRP SP SPIN 2.
 - c. Executive Series Centrifuge: 3.0 minutes and 4400 RPM.
 - d. Press the start button. Once the centrifuge stops, remove the **Device**.
39. SECOND EXTRACTION: Remove the clear cap from the bottom port. Using the 60mL syringe, aspirate plasma from the bottom port leaving 4mL or the desired amount in the device.
40. RESUSPEND THE PRP: Gently swirl the Concentrating Device to re-suspend the platelet concentrate into the plasma.
41. EXTRACT PRP: Attach a sterile 12mL syringe to the bottom port and tilt to immerse the aspirating pipe, then aspirate the platelet rich plasma through the open port of the aspirating pipe. Remove sterile syringe and apply a sterile cap.

INSTRUCTIONS FOR USE FOR 120mL SYSTEM













PREPARATION PROTOCOL:

42. NOTE: Use standard sterile aseptic technique throughout the following procedure. Always swab needle-less ports with alcohol before and after accessing.
43. WHOLE BLOOD DRAW: Attach the sterile filter needle onto two sterile 60mL syringes. Draw 8mL of Citrate Anticoagulant into each 60mL syringe. Remove the filter needle from the syringes. Attach the butterfly needle onto 60mL syringe and prime the needle with the anticoagulant. Slowly draw 52mL of whole blood per syringe, from the patient filling each syringe to 60mL. Gently, but thoroughly mix the blood and anticoagulant upon collection to prevent coagulation.

CONCENTRATING PROTOCOL:


42. LOAD: IMPORTANT: Attach sterile non-vented clear cap to the bottom port of EACH device. **STERILE CLEAR CAP MUST BE ATTACHED TO THE BOTTOM PORT BEFORE CENTRIFUGATION. WITH WHITE VENT OPEN**, slowly inject the anticoagulated whole blood through the top port of EACH Concentrating Device. **THEN CLOSE WHITE VENT ON EACH DEVICE.**
43. BALANCE: Make sure EACH device contains the same amount of volume. Then place them directly opposite to each other in the centrifuge rotor buckets. Close the lid.
44. FIRST SPIN:
 - a. Sapphire Series Centrifuge: PUREPRP 60 SPIN 1.
 - b. Platinum Series Centrifuge: PUREPRP SP SPIN 1.
 - c. Executive Series Centrifuge: 2.0 minutes and 4400 RPM.
 - d. Press the start button. Once the centrifuge stops, remove the **Device**.
45. FIRST EXTRACTION & TRANSFER: Attach the sterile 60mL syringe to the top port of EACH Device.
 - a. **LP-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate the remaining plasma buffycoat, leaving the RBC layer behind.
 - b. **MR-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate collecting 0.5mL of RBC buffycoat.
 - c. **NR-Protocol** - Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the AspiraDome™. Then, using a 3mL syringe, slowly aspirate collecting 1mL of RBC buffycoat..
 - d. **Remove CLEAR CAP** from bottom port and inject the 3mL syringe buffycoat solution, followed by the platelet plasma suspension (PPS). **Then SNUGGLY PLACE STERILE CLEAR CAP ON THE BOTTOM PORT!!!**
46. SECOND SPIN: Place Concentrating Devices back into the centrifuge rotor bucket directly opposite the counterbalance device. Close the lid.
 - a. Sapphire Series Centrifuge: PUREPRP 60 SPIN 2.
 - b. Platinum Series Centrifuge: PUREPRP SP SPIN 2.
 - c. Executive Series Centrifuge: 4.0 minutes and 4400 RPM.
 - d. Press the start button. Once the centrifuge stops, remove EACH **Device**.
47. SECOND EXTRACTION: Remove the clear cap from the bottom port of EACH device. Using a 60mL syringe, aspirate plasma from the bottom port, leaving 7mL or the desired amount in EACH device.
48. RESUSPEND THE PRP: Gently swirl the Concentrating Device to re-suspend the platelet concentrate into the plasma.
49. EXTRACT PRP: Attach a sterile 12mL syringe to the bottom port of EACH DEVICE and tilt to immerse the aspirating pipe, then aspirate the platelet rich plasma through the open port of the aspirating pipe. Remove sterile syringe and apply a sterile cap.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician

 Do not use if package is damaged	 Attention, read instruction for use	 Medical Device	 Single use only	 Store in a cool place	 Rx Only Prescription Use
 STERILEEO	 Do not re-sterilize	 Consult instruction for use	 Store in a dry place	 EmCyle Corporation 4331 Veronica S. Shoemaker Blvd. Fort Myers, FL 33916 Phone: 239-481-7725	 Non-pyrogenic

PREPARATION PROTOCOL


STEP 1:



Using the filtered needle, draw 8mL of Sodium Citrate Anticoagulant into 60mL Syringe. Then collect 52mL whole blood filling syringe to 60mL.

STEP 2:

**ATTACH
CLEAR
CAP**



Attach clear non vented cap to the bottom port.


CONCENTRATING PROTOCOL

STEP 3:



Inject anticoagulated whole blood through the top needle-less port.

STEP 4:



Counterbalance device with equal volume.



Place in the centrifuge rotor at opposite ends.


STEP 5:

**Sapphire Series Centrifuge:
PUREPRP 60 SPIN 1**


**Platinum Series Centrifuge:
PUREPRP SP SPIN 1**

**Executive Series Centrifuge
Set to 2 minutes and 4400**


STEP 6:




Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the Aspiradome™

LEUKOCYTE
POOR







Using a 3mL syringe, slowly aspirate the remaining plasma buffycoat, leaving the RBC layer behind.

MONOCYTE
RICH



Using a 3mL syringe, slowly aspirate collecting 0.5mL of RBC buffycoat. Fill until RBC layer reaches 0.5mL

NEUTROPHIL
RICH








Using a 3mL syringe, slowly aspirate collecting 1mL of RBC buffycoat. Fill until RBC layer reaches 1mL.







<p>STEP 7:</p>  <p>Inject the 3mL syringe buffycoat solution, followed by the platelet plasma suspension (PPS) through bottom port. Then place sterile clear cap.</p>	<p>STEP 8:</p> <p>Place back in centrifuge rotor. Process at:</p> <p>Sapphire Series Centrifuge: PUREPRP 60 SPIN 2</p> <p>Platinum Series Centrifuge: PUREPRP SP SPIN 2</p> <p>Executive Series Centrifuge 4 minutes and 4400</p>	<p>STEP 9:</p>  <p>After centrifugation, the platelet buffycoat will be separated at the bottom of the device.</p>	<p>STEP 10:</p>  <p>Connect the syringe to the bottom needle-less port and aspirate plasma, leaving 7mL or the desired amount in the device.</p>
<p>STEP 11:</p>  <p>Gently swirl to resuspend the platelet buffycoat into the plasma.</p>	<p>STEP 12:</p>  <p>Connect the 12mL syringe to the bottom needle-less port and tilt to aspirate the PRP</p>	<p>STEP 13:</p>  <p>7mL or desired amount PRP</p>	

PREPARATION PROTOCOL

<p>STEP 1:</p>  <p>Using the filtered needle, draw 3mL of Sodium Citrate Anticoagulant into 30mL syringe. Then collect 27mL whole blood filling syringe to 30mL.</p>	<p>STEP 2:</p>  <p>ATTACH CLEAR NON-VENTED CAP</p> <p>Attach clear non vented cap to the bottom port.</p>
--	--


CONCENTRATING PROTOCOL

<p>STEP 3:</p>  <p>Inject anticoagulated whole blood through the top needle-less port.</p>	<p>STEP 4:</p>  <p>Counterbalance device with equal volume.</p>  <p>Place in the centrifuge rotor at opposite ends.</p>	<p>STEP 5:</p> <p>Sapphire Series Centrifuge: PUREPRP 30 SPIN 1</p> <p>Platinum Series Centrifuge: PUREPRP SP SPIN 1</p> <p>Executive Series Centrifuge Set to 1 minutes and 4400</p>		
<p>STEP 6:</p>  <p>RBC layer at base of the Aspiradome™</p> <p>Aspirate the platelet plasma suspension (PPS) into the 30mL syringe, until the RBC interface reaches the bottom of the Aspiradome™</p>	<p>LEUKOCYTE POOR</p>  <p>Aspirate leaving RBC layer behind</p> <p>Using a 3mL syringe, slowly aspirate the remaining plasma buffycoat, leaving the RBC layer behind.</p>	<p>MONOCYTE RICH</p>  <p>Fill until RBC layer reaches 0.5mL</p> <p>Using a 3mL syringe, slowly aspirate collecting 0.5mL of RBC buffycoat. Fill until RBC layer reaches 0.5mL.</p>	<p>NEUTROPHIL RICH</p>  <p>Fill until RBC layer reaches 1mL</p> <p>Using a 3mL syringe, slowly aspirate collecting 1mL of RBC buffycoat. Fill until RBC layer reaches 1mL.</p>	

<p>STEP 7:</p>  <p>Inject the 3mL syringe buffycoat solution, followed by the platelet plasma suspension (PPS) through bottom port. Then place sterile clear cap.</p>	<p>STEP 8:</p> <p>Place back in centrifuge rotor. Process at:</p> <p>Sapphire Series Centrifuge: PUREPRP 30 SPIN 2</p> <p>Platinum Series Centrifuge: PUREPRP SP SPIN 2</p> <p>Executive Series Centrifuge 3 minutes and 4400</p>	<p>STEP 9:</p>  <p>After centrifugation, the platelet buffycoat will be separated at the bottom of the device.</p>	<p>STEP 10:</p>  <p>Leave 3-4mL or desired amount</p> <p>Connect the syringe to the bottom needle-less port and aspirate plasma, leaving 3-4mL or the desired amount in the device.</p>
<p>STEP 11:</p>  <p>Gently swirl to resuspend the platelet buffycoat into the plasma.</p>	<p>STEP 12:</p>  <p>Connect the 12mL syringe to the bottom needle-less port and tilt to aspirate the PRP</p>		<p>STEP 13:</p>  <p>3-4mL or desired amount PRP</p>


PREPARATION PROTOCOL

STEP 1:



Using the filtered needle, draw 8mL of Sodium Citrate Anticoagulant into two (2) 60mL Syringes. Then collect 52mL whole blood filling EACH syringe to 60mL.


STEP 2:



Attach clear non vented cap to the bottom port of EACH device.


CONCENTRATING PROTOCOL

STEP 3:




Inject anticoagulated whole blood into EACH device through the top needleless port.

STEP 4:



Counterbalance device with equal volume.



Place in the centrifuge rotor at opposite ends.


STEP 5:

Sapphire Series Centrifuge:
PUREPRP 60 SPIN 1

Platinum Series Centrifuge:
PUREPRP SP SPIN 1

Executive Series Centrifuge
Set to 2 minutes and 4400


STEP 6:



RBC layer at base of the Aspiradome™

FOR EACH DEVICE,
Aspirate the platelet plasma suspension (PPS) into the 60mL syringe, until the RBC interface reaches the bottom of the Aspiradome™


LEUKOCYTE POOR



Aspirate leaving RBC layer behind

Using a 3mL syringe, slowly aspirate the remaining plasma buffycoat, leaving the RBC layer behind.


MONOCYTE RICH



Fill until RBC layer reaches 0.5mL







Using a 3mL syringe, slowly aspirate collecting 0.05mL of RBC buffycoat. Fill until RBC layer reaches 0.05mL.

NEUTROPHIL RICH



Fill until RBC layer reaches 1mL

Using a 3mL syringe, slowly aspirate collecting 1mL of RBC buffycoat. Fill until RBC layer reaches 1mL.

<p>STEP 7:</p>  <p>Place sterile cap</p> <p>FORE EACH DEVICE. Inject the 3mL syringe buffycoat solution, followed by the platelet plasma suspension (PPS) through bottom port. Then place sterile clear cap.</p>	<p>STEP 8:</p> <p>Place back in centrifuge rotor at opposite ends.</p> <p>Process at:</p> <p>Sapphire Series Centrifuge: PUREPRP 60 SPIN 2</p> <p>Platinum Series Centrifuge: PUREPRP SP SPIN 2</p> <p>Executive Series Centrifuge 4 minutes and 4400</p>	<p>STEP 9:</p>  <p>Platelet buffycoat</p> <p>After centrifugation, the platelet buffycoat will be separated at the bottom of EACH device.</p>	<p>STEP 10:</p>  <p>Leave 7mL or desired amount</p> <p>Connect the syringe to the bottom needle-less port of EACH DEVICE and aspirate plasma, leaving 7mL or the desired amount in EACH device.</p>
<p>STEP 11:</p>  <p>Gently swirl EACH DEVICE to resuspend the platelet buffycoat into the plasma.</p>	<p>STEP 12:</p>  <p>Connect the 12mL syringe to the bottom needle-less port of EACH DEVICE and tilt to aspirate the PRP</p>	<p>STEP 13:</p>  <p>7mL or desired amount PRP</p>	